



Anti-Alpha 1-Acid Glycoprotein polyclonal antibody (CABT-BL8676)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Alpha 1-Acid Glycoprotein
Isotype	IgG
Source/Host	Goat
Species Reactivity	Canine, Dog
Conjugate	Unconjugated
Applications	ELISA, IP, WB
Format	Whole Serum
Size	1 ml
Preservative	None
Storage	4°C short term, freeze long term for 2 years