



Anti-Alpha 1-Microglobulin polyclonal antibody (CABT-BL8682)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Alpha 1-Microglobulin
Isotype	IgY
Source/Host	Chicken
Species Reactivity	Human
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Size	1 mg
Buffer	Phosphate buffered saline, pH 7.2
Preservative	0.1% Sodium Azide
Storage	4°C for 2 years