



Anti-Alpha 2-Macroglobulin polyclonal antibody (CABT-BL8684)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Alpha 2-Macroglobulin
Isotype	IgG
Source/Host	Goat
Species Reactivity	Rat
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, IP, WB
Format	Liquid
Size	100 µg
Buffer	Phosphate buffered saline, pH 7.2
Preservative	0.1% Sodium Azide
Storage	4°C for 2 years