



# Anti-Alpha Fetoprotein polyclonal antibody (CABT-BL8687)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                    |                                   |
|--------------------|-----------------------------------|
| Target             | Alpha Fetoprotein                 |
| Isotype            | IgG                               |
| Source/Host        | Goat                              |
| Species Reactivity | Human                             |
| Purification       | Antigen Affinity Purified         |
| Conjugate          | Unconjugated                      |
| Applications       | ELISA, WB                         |
| Format             | Liquid                            |
| Size               | 1 mg                              |
| Buffer             | Phosphate buffered saline, pH 7.2 |
| Preservative       | 0.1% Sodium Azide                 |
| Storage            | 4°C for 2 years                   |