



# Anti-beta Galactosidase (LacZ) polyclonal antibody (CABT-BL8692)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Target	beta Galactosidase
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, WB, IEP, IHC
Format	Liquid
Size	100 µg
Buffer	Phosphate buffered saline, pH 7.2
Preservative	0.1% Sodium Azide
Storage	4°C for 2 years