



Anti-C3 polyclonal antibody [Biotin] (CABT-BL8698)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	C3
Isotype	IgG
Source/Host	Goat
Species Reactivity	Guinea Pig
Purification	Antigen Affinity Purified
Conjugate	Biotin
Applications	ELISA, IHC, WB
Format	Liquid
Size	100 µg
Buffer	Phosphate buffered saline, pH 7.2 with 0.2% BSA
Preservative	0.1% Sodium Azide
Storage	4°C for 1 year

BACKGROUND

Introduction	Complement component C3 plays a central role in the activation of complement system. Its activation is required for both classical and alternative complement activation pathways. A peptide (C3a) derived from the encoded protein has antimicrobial activity, so people with C3 deficiency are susceptible to bacterial infection. [provided by RefSeq, Nov 2014]
--------------	---

GENE INFORMATION

Entrez Gene ID

[718](#)
