



Anti-Calmodulin Binding Peptide polyclonal antibody [HRP] (CABT-BL8712)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|--------------------|-----------------------------------|
| Target | Calmodulin Binding Peptide |
| Isotype | IgG |
| Source/Host | Rabbit |
| Species Reactivity | N/A |
| Purification | Antigen Affinity Purified |
| Conjugate | HRP |
| Applications | ELISA, IHC, WB |
| Format | Liquid |
| Size | 100 µg |
| Buffer | 50mM HEPES-Saline with stabilizer |
| Preservative | None |
| Storage | 4°C for 1 year |