



Anti-Canine Coronavirus polyclonal antibody (CABT-BL8713)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Canine Coronavirus
Isotype	IgG
Source/Host	Goat
Species Reactivity	Canine Coronavirus
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, IP, WB
Format	Liquid
Size	1 mg
Buffer	Phosphate buffered saline, pH 7.2
Preservative	0.1% Sodium Azide
Storage	4°C for 1 year