



Anti-Canine Heartworm polyclonal antibody [HRP] (CABT-BL8717)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Canine Heartworm
Isotype	IgY
Source/Host	Chicken
Species Reactivity	Canine Heartworm
Purification	Antigen Affinity Purified
Conjugate	HRP
Applications	ELISA, IHC, WB
Format	Liquid
Size	1 mg
Buffer	50mM HEPES-Saline with stabilizer
Preservative	None
Storage	4°C for 1 year