



Anti-CCDC75 polyclonal antibody (CABT-BL872)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|--------------------|---|
| Isotype | IgG |
| Source/Host | Rabbit |
| Species Reactivity | Human |
| Purification | This antibody is purified through a protein A column, followed by peptide affinity purification. |
| Conjugate | Unconjugated |
| Applications | WB, ELISA |
| Buffer | PBS with 0.09% (W/V) sodium azide |
| Preservative | 0.09% Sodium Azide |
| Storage | CCDC75 Antibody (N-term) can be refrigerated at 2-8°C for up to 6 months. For long term storage, place the at -2°C. |

GENE INFORMATION

| | |
|----------------|-----------------------------|
| Entrez Gene ID | 253635 |
| Protein Refseq | NP_777591.2 |
| UniProt ID | Q8N954 |