



# Anti-CCDC82 (aa 1-50) polyclonal antibody (CABT-BL874)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                    |                                                                         |
|--------------------|-------------------------------------------------------------------------|
| Immunogen          | Synthetic peptide from a region within amino acids 1-50 of Human CCDC82 |
| Isotype            | IgG                                                                     |
| Source/Host        | Rabbit                                                                  |
| Species Reactivity | Human                                                                   |
| Purification       | Immunogen affinity purified                                             |
| Conjugate          | Unconjugated                                                            |
| Applications       | WB, IP                                                                  |
| Format             | Liquid                                                                  |
| Buffer             | Tris citrate/phosphate, pH 7-8                                          |
| Preservative       | 0.09% Sodium Azide                                                      |
| Storage            | Store at +4°C.                                                          |

## GENE INFORMATION

|                |                            |
|----------------|----------------------------|
| Entrez Gene ID | <a href="#">79780</a>      |
| Protein Refseq | <a href="#">NP_079001</a>  |
| UniProt ID     | <a href="#">A0A024R3B3</a> |