



# Anti-CRP Neat polyclonal antibody (CABT-BL8744)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Target	CRP
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Dog
Conjugate	Unconjugated
Applications	ELISA, IP, WB
Format	Neat Antiserum
Size	1 ml
Preservative	0.1% Sodium Azide
Storage	4°C short term, freeze long term for 2 years

## BACKGROUND

### Introduction

The protein encoded by this gene belongs to the pentaxin family. It is involved in several host defense related functions based on its ability to recognize foreign pathogens and damaged cells of the host and to initiate their elimination by interacting with humoral and cellular effector systems in the blood. Consequently, the level of this protein in plasma increases greatly during acute phase response to tissue injury, infection, or other inflammatory stimuli. [provided by RefSeq, Sep 2009]

# GENE INFORMATION

Entrez Gene ID

[1401](#)

---