



# Anti-CRP polyclonal antibody [Biotin] (CABT-BL8746)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Target	CRP
Isotype	IgG
Source/Host	Goat
Species Reactivity	Human
Purification	Antigen Affinity Purified
Conjugate	Biotin
Applications	ELISA, IHC, WB
Format	Liquid
Size	1 mg
Buffer	Phosphate buffered saline, pH 7.2 with 0.2% BSA
Preservative	0.1% Sodium Azide
Storage	4°C for 1 year

## BACKGROUND

Introduction	The protein encoded by this gene belongs to the pentaxin family. It is involved in several host defense related functions based on its ability to recognize foreign pathogens and damaged cells of the host and to initiate their elimination by interacting with humoral and cellular effector systems in the blood. Consequently, the level of this protein in plasma increases greatly during
--------------	--

acute phase response to tissue injury, infection, or other inflammatory stimuli. [provided by RefSeq, Sep 2009]

---

## GENE INFORMATION

Entrez Gene ID [1401](#)

---