



Anti-6xHIS Tag polyclonal antibody [FITC] (CABT-BL8750)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	6xHIS tag
Isotype	IgY
Source/Host	Chicken
Species Reactivity	N/A
Purification	Antigen Affinity Purified
Conjugate	FITC
Applications	IF, IHC, WB
Format	Liquid
Size	100 µg
Buffer	Phosphate buffered saline, pH 7.2 with 0.2% BSA
Preservative	0.1% Sodium Azide
Storage	4°C for 1 year