



Anti-6xHIS Tag polyclonal antibody [HRP] (CABT-BL8751)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	6xHIS tag
Isotype	IgY
Source/Host	Chicken
Species Reactivity	N/A
Purification	Antigen Affinity Purified
Conjugate	HRP
Applications	ELISA, IHC, WB
Format	Liquid
Size	100 µg
Buffer	50mM HEPES-Saline with stabilizer
Preservative	None
Storage	4°C for 1 year