



Anti-HA Tag polyclonal antibody [FITC] (CABT-BL8763)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	HA tag
Isotype	IgY
Source/Host	Chicken
Species Reactivity	N/A
Purification	Antigen Affinity Purified
Conjugate	FITC
Applications	ELISA, IF, WB
Format	Liquid
Size	100 µg
Buffer	Phosphate buffered saline, pH 7.2 with 0.2% BSA
Preservative	0.1% Sodium Azide
Storage	4°C for 1 year