



# Anti-C-Myc-Agarose Immobilized polyclonal antibody (CABT-BL8773)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Target	C-Myc tag
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	IP: 0.1-1.0 mg
Format	Liquid
Size	200 µl
Buffer	Phosphate buffered saline, pH 7.2
Preservative	0.1% Sodium Azide
Storage	4°C for 1 year