



## Anti-DYKDDDDK polyclonal antibody [Biotin] (CABT-BL8783)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|                    |                                                 |
|--------------------|-------------------------------------------------|
| Target             | DYKDDDDK                                        |
| Isotype            | IgG                                             |
| Source/Host        | Rabbit                                          |
| Species Reactivity | N/A                                             |
| Purification       | Antigen Affinity Purified                       |
| Conjugate          | Biotin                                          |
| Applications       | ELISA, IHC, WB                                  |
| Format             | Liquid                                          |
| Size               | 100 µg                                          |
| Buffer             | Phosphate buffered saline, pH 7.2 with 0.2% BSA |
| Preservative       | 0.1% Sodium Azide                               |
| Storage            | 4°C for 1 year                                  |