



Anti-DYKDDDDK polyclonal antibody [HRP] (CABT-BL8787)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	DYKDDDDK
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Antigen Affinity Purified
Conjugate	HRP
Applications	ELISA, IHC, WB
Format	Liquid
Size	100 µg
Buffer	50mM HEPES-Saline with stabilizer
Preservative	None
Storage	4°C for 1 year