



# Anti-Elastase-IgG Fraction polyclonal antibody (CABT-BL8788)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Target	Elastase
Isotype	IgG
Source/Host	Sheep
Species Reactivity	Human
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Size	1 ml
Buffer	Phosphate buffered saline, pH 7.2
Preservative	0.1% Sodium Azide
Storage	4°C for 2 years