



# Anti-HA Tag polyclonal antibody [HRP] (CABT-BL8797)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Target	HA tag
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Antigen Affinity Purified
Conjugate	HRP
Applications	ELISA, IHC, WB
Format	Liquid
Size	100 µg
Buffer	50mM HEPES-Saline with stabilizer
Preservative	None
Storage	4°C for 1 year