



Anti-CCDC96 (aa 2-51) polyclonal antibody (CABT-BL880)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 2-51 (DVSSEHTKDP GGEGGDGSL AARPSKIKAS SGPPTSPEPG ELESEPEEEEE) of human CCDC96.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	257236
Protein Refseq	NP_699207
UniProt ID	Q2M329