



Anti-VSV-G Tag polyclonal antibody (CABT-BL8812)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	VSV-g tag
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, IEP, WB, IHC
Format	Liquid
Size	100 µg
Buffer	Phosphate buffered saline, pH 7.2
Preservative	0.1% Sodium Azide
Storage	4°C for 2 years