



Anti-CCDC99 (aa 458-488) polyclonal antibody (CABT-BL882)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide conjugated to KLH, corresponding to a region within C terminal amino acids 458-488 of Human CCDC99.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cytoplasm; cytoskeleton; centrosome. Chromosome; centromere; kinetochore. Nucleus. Cytoplasm; cytoskeleton; spindle pole. Localizes to the nucleus in interphase and to the kinetochore in early prometaphase. Relocalizes to the mitotic spindle pole be
Format	Liquid
Buffer	99% PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

GENE INFORMATION

Entrez Gene ID [54908](#)

Protein Refseq [NP_060255](#)

UniProt ID [Q96EA4](#)
