



Anti-GFP Antibody polyclonal antibody (CABT-BL8826)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Green Fluorescent Protein
Isotype	IgY
Source/Host	Chicken
Species Reactivity	N/A
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, IF, IHC, IP, WB
Format	Lyophilized
Size	0.5 mg
Buffer	Phosphate buffered saline, pH 7.2
Preservative	0.1% Sodium Azide
Storage	4°C for 2 years