



# Anti-CCHCR1 polyclonal antibody (CABT-BL884)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                              |   |
|------------------------------|---|
| <b>Immunogen</b>             | Synthetic peptide conjugated to KLH derived from within residues 750 to the C-terminus of Human CCHCR1.         |
| <b>Isotype</b>               | IgG   |
| <b>Source/Host</b>           | Rabbit  |
| <b>Species Reactivity</b>    | Human   |
| <b>Purification</b>          | Immunogen affinity purified   |
| <b>Conjugate</b>             | Unconjugated  |
| <b>Applications</b>          | WB  |
| <b>Cellular Localization</b> | Nuclear   |
| <b>Format</b>                | Liquid  |
| <b>Buffer</b>                | 1% BSA, PBS, pH 7.4   |
| <b>Preservative</b>          | 0.02% Sodium Azide  |
| <b>Storage</b>               | Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |

## BACKGROUND

|                     |  |
|---------------------|--|
| <b>Introduction</b> | CCHCR1 (Coiled coil alpha helical rod protein 1) is a regulator of keratinocyte proliferation or differentiation. It is overexpressed in keratinocytes in psoriatic lesions. |
|---------------------|--|

# GENE INFORMATION

Entrez Gene ID [54535](#)

---

Protein Refseq [NP\\_001099033](#)

---

UniProt ID [Q2TB68](#)

---