



Anti-CCHCR1 (center region) polyclonal antibody (CABT-BL885)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide (Human) selected from the center region, conjugated to KLH
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Cellular Localization	Nuclear
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

BACKGROUND

Introduction	CCHCR1 (Coiled coil alpha helical rod protein 1) is a regulator of keratinocyte proliferation or differentiation. It is overexpressed in keratinocytes in psoriatic lesions.
--------------	--

GENE INFORMATION

Entrez Gene ID [54535](#)

Protein Refseq [NP_001099033](#)

UniProt ID [Q2TB68](#)
