



Anti-IgY h+I polyclonal antibody [FITC] (CABT-BL8855)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	IgY h+I
Isotype	IgG
Source/Host	Goat
Species Reactivity	Chicken
Purification	Antigen Affinity Purified
Conjugate	FITC
Applications	IF, IHC, WB
Format	Liquid
Size	1 mg
Buffer	Phosphate buffered saline, pH 7.2 with 0.2% BSA
Preservative	0.1% Sodium Azide
Storage	4°C for 1 year