



Anti-IgY h+I Antibody polyclonal antibody (CABT-BL8856)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	IgY h+I
Isotype	IgG
Source/Host	Goat
Species Reactivity	Chicken
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, IP, LFIA(cap), WB
Format	Liquid
Size	1 mg
Buffer	Phosphate buffered saline, pH 7.2
Preservative	0.1% Sodium Azide
Storage	4°C for 2 years