



# Anti-IgY h+I Antibody polyclonal antibody [HRP] (CABT-BL8859)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Target	IgY h+I
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Chicken
Purification	Antigen Affinity Purified
Conjugate	HRP
Applications	ELISA, IHC, WB
Format	Liquid
Size	1 mg
Buffer	50mM HEPES-Saline with stabilizer
Preservative	None
Storage	4°C for 1 year