



# Anti-Kappa b+f polyclonal antibody [HRP] (CABT-BL8867)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Target	Kappa b+f
Isotype	IgG
Source/Host	Goat
Species Reactivity	Human
Purification	Antigen Affinity Purified
Conjugate	HRP
Applications	ELISA, IHC, WB
Format	Liquid
Size	1 mg
Buffer	50mM HEPES-Saline with stabilizer
Preservative	None
Storage	4°C for 1 year