



## Anti-NGAL polyclonal antibody (CABT-BL8882)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|                    |                                   |
|--------------------|-----------------------------------|
| Target             | NGAL (Lipocalin-2)                |
| Isotype            | IgG                               |
| Source/Host        | Goat                              |
| Species Reactivity | Rat                               |
| Purification       | Antigen Affinity Purified         |
| Conjugate          | Unconjugated                      |
| Applications       | ELISA, IP, WB                     |
| Format             | Liquid                            |
| Size               | 100 µg                            |
| Buffer             | Phosphate buffered saline, pH 7.2 |
| Preservative       | 0.1% Sodium Azide                 |
| Storage            | 4°C for 2 years                   |