



Anti-Plasmodium aldolase polyclonal antibody (CABT-BL8888)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Plasmodium aldolase
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Plasmodium
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, IP, WB
Format	Liquid
Size	1 mg
Buffer	Phosphate buffered saline, pH 7.2
Preservative	None
Storage	4°C for 2 years