



Anti-Transferrin polyclonal antibody (CABT-BL8910)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Transferrin
Isotype	IgG
Source/Host	Goat
Species Reactivity	Rabbit
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, IP, WB
Format	Liquid
Size	1 mg
Buffer	Phosphate buffered saline, pH 7.2
Preservative	0.1% Sodium Azide
Storage	4°C for 2 years