



Anti-CCV polyclonal antibody (CABT-BL8933)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Isotype	IgG
Source/Host	Pig
Species Reactivity	Canine Coronavirus
Conjugate	Unconjugated
Applications	IFA
Size	2 ml
Preservative	None
Storage	Store at 2-7 °C. Do not freeze!

BACKGROUND

Introduction	Canine Coronavirus is a single stranded RNA type of virus with a fatty protective coating. Because the virus is covered in a fatty membrane, it is relatively easily inactivated with detergent and solvent-type disinfectants. It is spread by virus shedding in the feces of infected dogs.
--------------	---