



Anti-CDV polyclonal antibody (CABT-BL8934)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Does not stain other common canine viruses.
Isotype	IgG
Source/Host	Goat
Species Reactivity	Canine Distemper Virus
Conjugate	Unconjugated
Applications	IFA
Size	2 ml
Preservative	None
Storage	Store at 2-7 °C. Do not freeze!

BACKGROUND

Introduction	Canine distemper is caused by a single-stranded RNA virus of the family Paramyxoviridae (the same family of the viruses that causes measles, mumps, and bronchiolitis in humans). The disease is highly contagious via inhalation and fatal 50% of the time. Despite extensive vaccination in many regions, it remains a major disease of dogs, and is the leading cause of infectious disease death in dogs.
--------------	---