



Anti-Clostridium CD1021 (N-terminal) polyclonal antibody (CABT-BL903)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to N terminal residues of Clostridium difficile protein CD1021
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	C. difficile
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Buffer	50% Glycerol, PBS
Preservative	0.01% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	4916761
UniProt ID	Q18AM0