



Anti-CD40LG polyclonal antibody (CABT-BL948)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide surrounding amino acid 246 of human CD40 Ligand
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity purified
Conjugate	Unconjugated
Applications	WB
Preservative	None
Storage	-20°C

BACKGROUND

Introduction	The protein encoded by this gene is expressed on the surface of T cells. It regulates B cell function by engaging CD40 on the B cell surface. A defect in this gene results in an inability to undergo immunoglobulin class switch and is associated with hyper-IgM syndrome. [provided by RefSeq, Jul 2008]
--------------	--

GENE INFORMATION

Entrez Gene ID	959
----------------	---------------------

Protein Refseq [NP_000065](#)

UniProt ID [P29965](#)
