



Anti-CD40LG polyclonal antibody [PE] (CABT-BL953)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Human IgM Fc5μ fragment
Source/Host	Donkey
Species Reactivity	Human
Purification	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Human IgM coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities, pepsin digestion and chromatographic separation. Ass
Conjugate	PE
Applications	IF
Format	Lyophilized
Buffer	0.01 M Sodium Phosphate, 0.25 M Sodium Chloride, pH 7.2
Preservative	None
Storage	2 years at -20°C

BACKGROUND

Introduction	The protein encoded by this gene is expressed on the surface of T cells. It regulates B cell function by engaging CD40 on the B cell surface. A defect in this gene results in an inability to undergo immunoglobulin class switch and is associated with hyper-IgM syndrome. [provided by RefSeq, Jul 2008]
---------------------	--

GENE INFORMATION

Entrez Gene ID [959](#)

Protein Refseq [NP_000065](#)

UniProt ID [P29965](#)
