



# Rabbit Anti-HPV 35-L1 Polyclonal antibody (CABT-CS038)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	This antibody has different degrees of cross-reactivity to L1 proteins of other types of HPV
<b>Target</b>	HPV 35 L1
<b>Immunogen</b>	HPV 35 L1
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	HPV
<b>Purification</b>	Protein G purified
<b>Conjugate</b>	unconjugated
<b>Applications</b>	WB, ELISA
<b>Format</b>	Liquid ascites fluid
<b>Concentration</b>	1 mg/mL
<b>Size</b>	100 µg
<b>Buffer</b>	Tris-Gly
<b>Preservative</b>	None
<b>Storage</b>	Store at -20°C. Avoid repeated freezing and thawing cycles.
<b>Ship</b>	Shipped with wet ice

# BACKGROUND

Introduction	Human papillomavirus infection (HPV infection) is an infection caused by the human papillomavirus (HPV). About 90% of HPV infections cause no symptoms and resolve spontaneously within two years. However, in some cases, an HPV infection persists and results in either warts or precancerous lesions.
Keywords	HPV 35 L1; HPV; HPV 35; HPV L1