



Rabbit Anti-VSV Polyclonal antibody (CABT-CS059)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	VSV
Target	VSV
Immunogen	Purified stock of VSV-hIFN β virus
Isotype	Serum
Source/Host	Rabbit
Species Reactivity	VSV
Conjugate	unconjugated
Applications	WB, IF, IHC, FC, ELISA, Neut
Format	Liquid
Size	0.25 ml
Preservative	None
Storage	Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

BACKGROUND

Introduction	Vesicular stomatitis virus (VSV) – A virus in the family Rhabdoviridae (like Rabies), genus Vesiculovirus. Of the numerous vesiculoviruses, at least six have been confirmed to cause animal or human disease. It has particular importance to farmers in certain regions of the world where it infects cattle. This is because its clinical presentation is identical to the very important
---------------------	--

foot and mouth disease virus.

Keywords

Vesicular stomatitis virus; VSV; VSIV
