



Human Anti-SARS-CoV-2 Nucleoprotein Monoclonal antibody (purified), clone CON12 (CABT-CS061)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Binding to N-ter of SARS-CoV-2 NP
Target	SARS-CoV-2 NP
Isotype	IgM
Source/Host	Humanized
Species Reactivity	SARS-CoV-2
Clone	CON12
Purification	Anti-IgM column. ≥85% as determined by SDS-PAGE
Conjugate	unconjugated
Applications	ELISA
Format	Liquid
Concentration	1 mg/mL(may vary lot by lot)
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

BACKGROUND

Introduction

SARS-CoV-2 is the virus responsible for COVID-19. It was previously known as 2019 Novel coronavirus (2019-nCoV), is a positive-sense single-stranded RNA virus. SARS-CoV-2 has close genetic similarity to bat coronaviruses, from which it likely originated. An intermediate animal reservoir such as a pangolin is also thought to be involved in its introduction to humans. From a taxonomic perspective, SARS-CoV-2 is classified as a strain of the species severe acute respiratory syndrome-related coronavirus (SARSr-CoV).

Keywords

SARS-CoV-2; coronavirus; SARS-CoV-2 NP; SARS-CoV-2 Nucleocapsid Protein