



Rabbit Anti-MDMV Polyclonal antibody (CABT-CS165)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	There is no cross reaction with healthy plant tissues such as corn and wheat.
Target	MDMV
Immunogen	MDMV
Source/Host	Rabbit
Species Reactivity	MDMV
Conjugate	unconjugated
Applications	ELISA
Format	Liquid
Size	500 tests
Buffer	PBS
Storage	Store at 4°C.

BACKGROUND

Keywords	MDMV; Maize Dwarf Mosaic Virus
----------	--------------------------------