



Rabbit Anti-MSV Polyclonal antibody (CABT-CS166)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	There is no cross reaction with healthy plant tissues such as corn and wheat.
Target	MSV
Immunogen	MSV
Source/Host	Rabbit
Species Reactivity	MSV
Conjugate	unconjugated
Applications	ELISA
Format	Liquid
Size	500 tests
Buffer	PBS
Storage	Store at 4°C.

BACKGROUND

Keywords	MSV; Maize Streak Virus
----------	-------------------------