



# Rabbit Anti-YMV Polyclonal antibody (CABT-CS167)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Specificity	There is no cross reaction with healthy plant tissues such as rice and other host plants.
Target	YMV
Immunogen	YMV
Source/Host	Rabbit
Species Reactivity	YMV
Conjugate	unconjugated
Applications	ELISA
Format	Liquid
Size	500 tests
Buffer	PBS
Storage	Store at 4°C.

## BACKGROUND

Keywords	YMV; Yam Mosaic Virus
----------	-----------------------