



Rabbit Anti-PSV Polyclonal antibody (CABT-CS170)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	There is no cross reaction with healthy plant tissues such as peanut, bean, pea, tobacco, or other host plants.
Target	PSV
Immunogen	PSV
Source/Host	Rabbit
Species Reactivity	PSV
Conjugate	unconjugated
Applications	ELISA
Format	Liquid
Size	500 tests
Buffer	PBS
Storage	Store at 4°C.

BACKGROUND

Keywords PSV; Peanut Stunt Virus