



Rabbit Anti-POTY Polyclonal antibody (CABT-CS171)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	There is no cross reaction with healthy plant tissues.
Target	POTY
Immunogen	POTY
Source/Host	Rabbit
Species Reactivity	POTY
Conjugate	unconjugated
Applications	ELISA
Format	Liquid
Size	500 tests
Buffer	PBS
Storage	Store at 4°C.

BACKGROUND

Keywords	POTY; Potyvirus Group; Potyvirus
----------	----------------------------------