



Mouse Anti-HPV6 L1 Monoclonal antibody, clone 2642 (CABT-CS174)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Reacts with the major capsid protein L1
Target	HPV6 L1
Immunogen	HPV6 antibody was raised in mouse using papilloma virus type 6 as the immunogen.
Source/Host	Mouse
Species Reactivity	HPV6
Clone	2642
Conjugate	unconjugated
Applications	ELISA
Format	Liquid
Size	200 µg
Buffer	0.01M PBS, pH 7.3
Preservative	0.1% Sodium Azide
Storage	Aliquot and store at -20°C. Avoid repeated Freeze/Thaw cycles.

BACKGROUND

Introduction	Infection with specific types of human Papillomavirus (HPV) has been associated with an increased risk of developing cervical neoplasia. HPV types 6 and 11 have been associated with
--------------	---

relatively benign diseases such as genital warts but types 16 and 18 are the causative agents of cervical, vaginal and vulvar malignancies.

Keywords

Human papillomavirus type 6; HPV6; HPV 6; HPV6 L1; HPV 6 L1
