



Mouse Anti-S. Suis 2 Polyclonal antibody (CABT-CS178)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Streptococcus Suis 2
Target	Streptococcus Suis 2
Immunogen	whole cell protein of Streptococcus suis Type 2
Isotype	IgG
Source/Host	Mouse
Species Reactivity	S. Suis
Conjugate	unconjugated
Applications	ELISA, IF, IHC-P, WB
Format	Liquid
Concentration	1 mg/ml
Size	100 µg, 200 µg
Buffer	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Preservative	None
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

BACKGROUND

Keywords

Streptococcus Suis 2; Streptococcus Suis; Streptococcus; S. Suis
