



# Rabbit Anti-H7N9 HA Polyclonal antibody (CABT-CS179)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	Has cross-reactivity in ELISA with H7N9 (A/Anhui/1/2013) HA
<b>Target</b>	H7N9 HA
<b>Immunogen</b>	Recombinant Influenza A H7N9 (A/Shanghai/1/2013) HA protein.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	H7N9
<b>Purification</b>	Protein A & Antigen Affinity
<b>Conjugate</b>	unconjugated
<b>Applications</b>	ELISA, WB
<b>Format</b>	Liquid
<b>Size</b>	50 µl, 100 µl
<b>Buffer</b>	PBS
<b>Preservative</b>	None
<b>Storage</b>	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

## BACKGROUND

**Keywords**

nfluenzavirus A; Influenza A virus; Influenza A virus H7N9 HA; H7N9 HA; IAV H7N9 HA; H7N9; H7N9 Hemagglutinin

---